


# Multiplexing

---



# Pooling Libraries from Different Samples

Index sequences are “barcodes” for multiplexing




- Multiplexing involves pooling libraries from different biological samples to be sequenced together on the same flow cell.
- The i5 and i7 index sequences are barcodes that are shared by all molecules from the same library so that libraries can be distinguished from each other during data analysis.

# General Multiplexing Recommendations

- Incorporate library-specific barcodes. Avoid home-brew methods that add the barcode directly at the end of the DNA insert unless you know what you are doing.



- Be conservative about pooling
  - Not all libraries will be equally represented in your mix
  - Number of clusters may be lower than anticipated
- Consider how cross-contamination will affect your analysis and use redundant dual-indexing when possible.



# Single vs Dual Indexing

Figure 1 Single-Indexed Sequencing





Figure 2 Dual-Indexed Single-Read Sequencing



Dual indexing can either reduce library cross-contamination if indexes are used in a redundant fashion...  
...or increase the degree of multiplexing if indexes are used in a combinatorial fashion.

# Index “Hopping” and Library Cross-Contamination

Even with perfect lab technique, library contamination occurs on the flow cell.





Table 1: Best Practices for Reducing Index Hopping

Mitigation/Recommendation	Benefit/Outcome
Prepare dual indexed libraries with unique indexes <sup>a</sup>	Converts index hopped reads to undetermined
Sequence one 30× human genome per lane <sup>b</sup>	Avoids pooling and index hopping
Remove adapters (cleanup, spin columns, etc) <sup>c</sup>	Reduces levels of index hopping
Store prepared libraries at recommended temperature of −20° C°	Reduces levels of index hopping
Pool similar RNA-Seq samples together	Reduces contamination between high and low-expressors

Image Sources: Illumina; Sinha et al. 2017 bioRxiv